



CRISPR-Cas12a DNA detection kit (2-step)(lyophilized)(paper strip)

✉ info@ezassay.com

🌐 www.ezassay.com

EZassay Biotechnology Ltd.


Catalog Code:D-L-CAS12-LYO-2S

CONTENTS

<u>Contents</u>	<u>Page</u>
Brief introduction	1
Materials supplied	1
Other materials required	1
Storage	2
Sample for detection	2
Assay procedure	2
Interpretation of results	4
Schematic Diagram of Principle	5
Notes	6

Brief introduction

This kit combines isothermal amplification and CRISPR/Cas12a to achieve sensitive and specific DNA detection. This kit contains the universal reagents for both the isothermal amplification and the CRISPR/Cas12a detection. It makes R&D of point-of-care products easily.



Materials supplied

#	Component	Volume
1	Reaction Tube	96T
2	Reaction Buffer (2X)	1000 μ l
3	Starter (10X)	200 μ l
4	Cleavage Buffer (10X)	240 μ l
5	Diluent	8000 μ l
6	Cas12a Protein_ lyophilized	Add PCR water according to the label on the tube before use.
7	Cas12a Reporter_ lyophilized	
8	Positive Control_ lyophilized (primer and DNA template included)	
9	crRNA for Positive Control_ lyophilized	

#	Component	Volume
01	Reaction Buffer (2X)	10 μ l
02	Forward Primer (20 μ M) Reverse Primer (20 μ M)	0.5 μ l 0.5 μ l
03	DNA template*	x μ l
04	Starter (10X) **	2 μ l
05	Nuclease-free H ₂ O	To a total volume of 20 μ l

*For no template control group, use Nuclease-free water. For positive control group, add 2 μ L Positive Control as provided in the kit.

** Add starter at the final step and mix thoroughly.

- Flick to mix, quick spin down. Repeat 3 times.
- Incubate at 39°C for 20~40 minutes.
- Run agarose gel to make sure the amplification works.

2. The second part is CRISPR/Cas12a detection. Set the working temperature at 37 °C.

- Take 20 μ l reaction volume as an example. Prepare reactions as described in the table below. Suggest set up the reaction on ice.

#	Component	Volume
01	Cleavage Buffer (10X)	2 μ l
02	Cas12a Reporter (4 μ M)	0.25 μ l
03	Cas12a Protein (1 μ M) *	1 μ l
04	crRNA (Cas12a) (1 μ M) **	1 μ l
05	Amplicons from last step ***	x μ l
06	Nuclease-free H ₂ O	To a total volume of 20 μ l

* Diluent Cas12a Protein properly with ddH2O before use.

**The positive control group adds 1 μ L of crRNA for Positive Control (20X). For the other group, add the specific crRNA for sequence of interest.


*** $x \leq 5\mu$ L.

- Flick to mix, quick spin down. Repeat 3 times. Incubate tube at 37 °C for 30 ~ 60 minutes.

3. Observe the result.

- Take 1~10 μ L of the product and dilute it 9~200 folds with Diluent. (e.g. Add 2uL product into 78uL Diluent.)
- Take 70 μ L of the dilution and add it to paper strip, and observe the result within 5 minutes.


Interpretation of results



Lateral flow paper strip (indirect method for CRISPR SHERLOCK)

Catalog number: CS-FMBO-96.

Schematic Diagram of Principle



Attention

- l Designate and use distinct areas for sample preparation, reaction setup, and analysis to avoid carry-over contamination.
- l Do not open the tube after amplification if possible. Please enquiry EZassay team for technical support.
- l If PCR thermal cycler is used, please make sure to turn off the heat lid function or set to 40 °C.
- l The optimal concentration of primer, Cas protein, and crRNA varies depending on the projects. Optimized range of working concentration: primers: from 300nM~800nM for each primer; crRNA: 20nM~1000nM; reporter: 20nM~1000nM; Cas12a: 20nM~200nM.