



# RT-LAMP-CRISPR RNA detection kit (one-step)

✉ [info@ezassay.com](mailto:info@ezassay.com)

🌐 [www.ezassay.com](http://www.ezassay.com)

EZassay Biotechnology Ltd.

Catalog Code: RT-LAM-CAS-01  
RT-LAM-CAS-10

# CONTENTS

| <u>Contents</u>         | <u>Page</u> |
|-------------------------|-------------|
| Product Introduction    | 1           |
| Materials supplied      | 1           |
| Required materials      | 1           |
| Storage                 | 1           |
| Heat Inactivation       | 1           |
| Preparation of reagents | 3           |
| Notes                   | 3           |

## Product Introduction


This product combines RT-LAMP and CRISPR/Cas12b to achieve sensitive and specific RNA detection. All reactions are performed in one tube and in one step. It makes R&D of point-of-care products easily.

## Materials supplied

| Component                   | RT-LAM-CAS-01<br>96T | RT-LAM-CAS-10<br>96T*10 |
|-----------------------------|----------------------|-------------------------|
| Reaction Buffer (2X)        | 1.5ml                | 15 ml                   |
| RT-LAMP Enzyme Mix (25X)    | 100 $\mu$ l          | 1000 $\mu$ l            |
| Cas12b Protein (10 $\mu$ M) | 80 $\mu$ L           | 800 $\mu$ l             |
| Reporter (20 $\mu$ M)       | 130 $\mu$ l          | 1300 $\mu$ l            |
| Positive Control            | 50 $\mu$ l           | 500 $\mu$ l             |

## Required materials

- Template RNA
- LAMP primers (online software: <http://primerexplorer.jp/e/>)
- Q-PCR machine or heat block
- Nuclease-free water
- sgRNA for AapCas12b. Recommend High Yield crRNA/sgRNA Synthesis Kit (Cat.: SG-RNA-001)  
(AapCas12b crRNA scaffold sequence: 5' -  
GUCUAGAGGACAGAAUUUUUCAACGGGUGUGCCAAUGGCCACUUUCCAGGUGGCAAAGCCCG  
UUGAGCUUCUCAAUCUGAGAAGUGGCAC -3')



## Storage

-20°C. Suggest aliquot. Avoid repeated freeze-thaw.

## Preparation of reagents

1. Mix the reagents after melting on ice.
2. Set Q-PCR working temperature at 60 °C. Turn off the lid heating function or set to 60 °C.
3. Prepare reactions as described in the table below.


| Component                | Volume     | Working concentration                             |
|--------------------------|------------|---|
| Reaction Buffer (2X)     | 12.5µl     | 1X  |
| RT-LAMP Enzyme Mix (25X) | 1µl        | 1X  |
| Cas12b Protein (10µM)    | 0.75µL     | 300nM   |
| Primers (10X)            | 1.25µL     | 1.6 µM FIP/BIP,<br>0.4 µM LF/LB,<br>0.2 µM F3/B3, |
| Reporter (20µM)          | 1.25µL     | 1µM   |
| sgRNA (25µM)             | 1µL        | 1µM   |
| Template DNA*            | 1µL        | >10 copies or more                                |
| ddH2O                    | Up to 25µL | -   |

- \* For no template control group, use Nuclease-free water.  
For positive control group, add 5µL Positive control that provided in the kit.  
Template, primers, and probe with FAM label are already mixed as positive control.
4. Read fluorescence in Q-PCR machine for 30 to 60 minutes at 60°C.

## Notes

- Designate and use distinct areas for sample preparation, reaction setup, and analysis to avoid carry-over contamination.
- It is recommended to add a negative control without template for each experiment.
- To guarantee good reproducibility, it is recommended that the template RNA be added last.
- If Q-PCR machine is used, please make sure to turn off the heat lid function or set to 60 °C.
- RNase inhibitor (Cat. MRI-2000) and RNase-free consumables are recommended for RNA experiment.
- The optimal concentration of sgRNA and reporter varies depending on the projects.

## Nucleic acid detection based on CRISPR/Cas technology diagram



© Ezassay Biotechnology, Inc 03

Professional supplier of point-of-care test products

**EZ**assay **EZassay Biotechnology Ltd.**

[www.ezassay.com](http://www.ezassay.com)

[info@ezassay.com](mailto:info@ezassay.com)