



**Lateral flow paper strip for  
nucleic acid detection  
Instructions(for RPA)**

✉ [info@ezassay.com](mailto:info@ezassay.com)

🌐 [www.ezassay.com](http://www.ezassay.com)

EZassay Biotechnology Ltd.


Catalog Code: PS-FMBO-48  
PS-FMBO-96

# CONTENTS

<u>Contents</u>	<u>Page</u>
Product Introduction	1
Materials supplied	1
Procedure	1
Result	2
Note	2

## Product Introduction

This product detects FITC/FAM and biotin dual-labeled gene amplification products based on lateral flow technology using gold particles. One primer/probe labeled with biotin and the other with fluorescein isothiocyanate (FITC) or 6-carboxyfluorescein (6-FAM), ensuring that both labels can be incorporated into the double-stranded amplification product at the same time. Recommend to use with EZassay™ products Isothermal DNA amplification Nfo kit (Cat.#: NF-LYO-96)



## Materials supplied

Item	PS-FMBO-48	PS-FMBO-96	Store
Lateral flow paper strip (without cassette)	48 tests	96 tests	Room temperature (15- 30°C). Do not freeze. Protected from moisture and direct sunlight.

## Procedure

1. Diluent amplification product as needed using diluent buffer or ddH<sub>2</sub>O. (Experiments are required to determine the optimal dilution factor. Normally diluent 10~20 times for RPA amplification product. In some tests, for example, LAMP, no dilution is needed.)
2. Transfer 70µL diluent product to sample pad with a pipette. Or simply insert strip into the liquid. (The liquid level should not exceed the MAX line.)
3. Read result in 5 minutes.

## Result




Figure 1: Schematic diagram of the interpretation of the results

## Notes

- 1) Pay attention to carry-over contamination. Suggest performing amplification and observation in different areas. Contact EZassay biotech. technician for solutions of carry-over contamination.
- 2) Keep sealed. Avoid moisture.
- 4) Please put the used test strip in a sealed bag and dispose of it properly.